

**Aumictech**

Test & measurement • bought, sold, repaired, calibrated

Data Sheet

August 13, 2026

STANFORD RESEARCH SYSTEMS

Stanford Research QCM100

Stanford Research Systems QCM100 Quartz Crystal Microbalance Analog Controller The SRS QCM100 is an analog quartz crystal microbalance controller used with the QCM25 5 MHz crystal oscillator. It measures mass and viscoelastic loading by providing a 5 MHz TTL frequency output and analogarithmic conductance voltage related to series resistance. Unlike the QCM200 digital controller, the QCM100 requires an external frequency counter and precision voltmeter. Automatic gain control sustains series resonance under heavy loads up to 5000 ohm, with shunt-capacitance cancellation from 10 pF to 40 pF. The SRS QCM100 is an analog quartz crystal microbalance controller used with the QCM25 5 MHz crystal oscillator. It measures mass and viscoelastic loading by providing a 5 MHz TTL frequency output and analogarithmic conductance voltage related to series resistance. Unlike the QCM200 digital controller, the QCM100 requires an external frequency counter and precision voltmeter. Automatic gain control sustains series resonance under heavy loads up to 5000 ohm, with shunt-capacitance cancellation from 10 pF to 40 pF. Thin-film adsorption and desorption; electrochemical QCM (EQCM) with a potentiostat via crystal-face bias; liquid and polymer viscosity/elasticity studies; real-time nanogram-class mass monitoring; teaching and OEM setups that already have a counter and DVM.



MODEL

Stanford Research Systems QCM100

CALIBRATION

Available on request

CONDITION

Used

STATUS

In stock

RFQ / SKU

QCM100

LISTING

<https://aumictechlabs.com/products/qcm100-2/>

NIST-traceable cal

SAM registered

30-day returns

NAICS 423690 • 811210 • 541380

This unit

The Stanford Research Systems QCM100 is a quartz crystal microbalance (QCM) instrument that measures mass and viscosity changes at surfaces and within thin films in real time. Operating at 5 MHz with an AT-cut quartz resonator, it detects mass variations from fractions of a nanogram to micrograms, with sensitivity extending to sub-monolayers of atoms. The system quantifies both frequency shifts - proportional to deposited mass - and changes in series resonance resistance, which correlate directly to material viscosity and elasticity. A Q-factor exceeding 60,000 in liquid environments ensures measurement stability even in demanding chemical and biological applications.

Technical Specifications

Resonator Frequency: 5 MHz nominal, AT-cut plano-plano geometry
Crystal Diameter: 1 inch
Mass Detection Range: Fractions of a nanogram to micrograms
Detection Limit: Sub-monolayers of atoms
Heavy Load Capability: Maintains oscillation up to 5 k Ω resistance
Aqueous Stability: Sustained operation in solutions exceeding 88% glycerol (w/w)
Q-Factor: >60,000 typical in liquid environments
Crystal Surface Finish: Optical polish (<10 Å RMS)
Electrode Materials: Cr/Au standard; Ti/Au, Ti/Pt, and ITO optional
Frequency Output: 2.4 Vpp square wave, 50 Ω , TTL levels
Conductance Output (Vc): 0 to 10.625 VDC, 1 k Ω source impedance

Key Features

Measures resonant frequency shifts linearly proportional to surface mass deposition
Converts resistance changes to viscosity and elasticity data via: $R = 10,000 \times (10 - V_c/5) - 75 \Omega$
Flat (non-planoconvex) crystal geometry for uniform measurement fields
Compatible with liquid and gas environments through compact crystal holder
Supports alternative frequencies (6 MHz, 9 MHz) and specialized electrode coatings

Typical Applications

Thin-film deposition monitoring
Corrosion studies
Biosensing and biofilm characterization
Chemical sensing and detection

Compatibility & Integration

The complete system comprises the QCM100 analog controller, QCM25 crystal oscillator, crystal holder, and three 5 MHz quartz crystals. External instrumentation required: precision frequency counter for frequency measurement and high-resolution digital voltmeter (minimum 6-digit resolution with computer interface) for conductance measurement.

Features & description

From manufacturer datasheet

Features

- Analog controller plus QCM25 oscillator, Kynar holder, and three 5 MHz Cr/Au crystals
- Frequency output isaTTL square wave (about 2.4 Vpp into 50 ohm) for any commercial counter
- Conductance voltage V_c maps series resistance: $R = 10000 \times (10^{(-V_c/5)}) - 75 \text{ ohm}$
- Ten-turn capacitance-cancellation dial and Adjust-mode lock-in null of C_o
- Same 5 kOhm heavy-load capability as QCM200
- Cat-5 cable (about 3 ft) between controller and QCM25

QCM100 Characteristics / Specifications

PARAMETER	VALUE
Model	QCM100
Manufacturer	Stanford Research Systems
Instrument Type	Quartz Crystal Microbalance Analog Controller
Companion Oscillator	QCM25 5 MHz crystal oscillator
Alias	QCM100; SRS QCM100
Crystal Frequency	5 MHz, AT-cut, plano-plano
Crystal Diameter	1 inch
Electrode Materials	chromium/gold standard (polished/unpolished); Ti/Au, Ti/Pt, indium tin oxide optional
Frequency Output Frequency	5 MHz nominal
Frequency Output Level	TTL square wave; about 2.4 Vpp into 50 ohm
Frequency Output Impedance	50 ohm
Frequency Output Connector	BNC
Resistance Calculation	$R = 10000 \times (10^{(-V_c/5)}) - 75$ ohm
Resistance Range	0 to 5000 ohm
Conductance Voltage Level	0 to 10.625 VDC, logarithmic scale
Conductance Output Impedance	1 kOhm
Conductance Output Connector	BNC
Capacitance Cancellation Range	10 pF to 40 pF (20 pF nominal)
Achievable Cancellation	0.01 pF
Typical Parasitic Capacitance Co	about 20 pF (electronics on holder; no cable C)
Typical Motional Capacitance Cm	33 fF (1 inch 5 MHz crystal)
Typical Motional Inductance Lm	30 mH
Typical Motional Resistance Dry	10 ohm
Typical Motional Resistance Water	400 ohm (one face in water)
Typical Motional Resistance Glycerol	3500 ohm (one face in 88 % glycerol)

PARAMETER	VALUE
Sauerbrey Sensitivity Factor	56.6 Hz per ug per cm^2 (5 MHz AT-cut at room temperature)
Frequency Mass Coefficient	0.057 Hz per (ng/cm^2) (theory note)
Crystal Holder Material	KynarORing Material: Viton
Holder Connector	BNC
Controller Oscillator Cable	Cat-5, approximately 3 ft
Line Voltage	90 to 132 VAC or 175 to 264 VAC, 47 to 63 Hz
Power Consumption	20 VA total
Operating Mode	analog; external counter and DVM required (not included)
Heavy Load Capability	up to 5000 ohm series resistance
Dry Dial Preset	ten-turn dial 8.0 (dry operation setting per setup procedure)
Varactor Fixed Bias Preset	ten-turn dial 5.0 yields about 6.0 V reverse bias, Cv about 18 pF
Adjust Mode Modulation	75 Hz sine on varactor bias with High/Low/Nullled LEDs
Standard Crystals Included	three chromium/gold QCM crystals Notes Specs from SRS QCM100/QCM25 operation manual, QCM100/QCM200 combined catalog, and QCM theory/calibration note. Do not invent counter resolution: frequency accuracy is that of the user-supplied counter and timebase. Forabuilt-in 0.01 Hz counter use QCM200.
Conductance voltage Vc maps series resistance	$R = 10000 \times (10^{(-Vc/5)}) - 75 \text{ ohm}$

Full OEM specifications

Complete technical content extracted from the manufacturer datasheet. Sparse OEM tables are reconstructed as tables below; other text is preserved as paragraphs.

Specifications	
Technical Specifications	
PARAMETER	VALUE
Model	QCM100
Manufacturer	Stanford Research Systems
Instrument Type	Quartz Crystal Microbalance Analog Controller
Companion Oscillator	QCM25 5 MHz crystal oscillator
Alias	QCM100; SRS QCM100
Crystal Frequency	5 MHz, AT-cut, plano-plano
Crystal Diameter	1 inch
Electrode Materials	chromium/gold standard (polished/unpolished); Ti/Au, Ti/Pt, indium tin oxide optional
Frequency Output Frequency	5 MHz nominal
Frequency Output Level	TTL square wave; about 2.4 Vpp into 50 ohm
Frequency Output Impedance	50 ohm
Frequency Output Connector	BNC
Resistance Calculation	$R = 10000 \times (10^{(-Vc/5)}) - 75 \text{ ohm}$
Resistance Range	0 to 5000 ohm
Conductance Voltage Level	0 to 10.625 VDC, logarithmic scale
Conductance Output Impedance	1 kOhm
Conductance Output Connector	BNC
Capacitance Cancellation Range	10 pF to 40 pF (20 pF nominal)
Achievable Cancellation	0.01 pF
Typical Parasitic Capacitance Co	about 20 pF (electronics on holder; no cable C)
Typical Motional Capacitance Cm	33 fF (1 inch 5 MHz crystal)
Typical Motional Inductance Lm	30 mH
Typical Motional Resistance Dry	10 ohm
Typical Motional Resistance Water	400 ohm (one face in water)
Typical Motional Resistance Glycerol	3500 ohm (one face in 88 % glycerol)

PARAMETER	VALUE
Sauerbrey Sensitivity Factor	56.6 Hz per ug per cm^2 (5 MHz AT-cut at room temperature)
Frequency Mass Coefficient	0.057 Hz per (ng/cm^2) (theory note)
Crystal Holder Material	Kynar
O Ring Material	Viton
Holder Connector	BNC
Controller Oscillator Cable	Cat-5, approximately 3 ft
Line Voltage	90 to 132 VAC or 175 to 264 VAC, 47 to 63 Hz
Power Consumption	20 VA total
Operating Mode	analog; external counter and DVM required (not included)
Heavy Load Capability	up to 5000 ohm series resistance
Dry Dial Preset	ten-turn dial 8.0 (dry operation setting per setup procedure)
Varactor Fixed Bias Preset	ten-turn dial 5.0 yields about 6.0 V reverse bias, Cv about 18 pF
Adjust Mode Modulation	75 Hz sine on varactor bias with High/Low/Nullled LEDs
Standard Crystals Included	three chromium/gold QCM crystals

Specifications

PARAMETER	VALUE
Model	QCM100
Manufacturer	Stanford Research Systems
Instrument Type	Quartz Crystal Microbalance Analog Controller
Companion Oscillator	QCM25 5 MHz crystal oscillator
Alias	QCM100; SRS QCM100
Crystal Frequency	5 MHz, AT-cut, plano-plano
Crystal Diameter	1 inch
Electrode Materials	chromium/gold standard (polished/unpolished); Ti/Au, Ti/Pt, indium tin oxide optional
Frequency Output Frequency	5 MHz nominal
Frequency Output Level	TTL square wave; about 2.4 Vpp into 50 ohm
Frequency Output Impedance	50 ohm
Frequency Output Connector	BNC
Resistance Calculation	$R = 10000 \times (10^{(-Vc/5)}) - 75$ ohm
Resistance Range	0 to 5000 ohm
Conductance Voltage Level	0 to 10.625 VDC, logarithmic scale

PARAMETER	VALUE
Conductance Output Impedance	1 kOhm
Conductance Output Connector	BNC
Capacitance Cancellation Range	10 pF to 40 pF (20 pF nominal)
Achievable Cancellation	0.01 pF
Typical Parasitic Capacitance Co	about 20 pF (electronics on holder; no cable C)
Typical Motional Capacitance Cm	33 fF (1 inch 5 MHz crystal)
Typical Motional Inductance Lm	30 mH
Typical Motional Resistance Dry	10 ohm
Typical Motional Resistance Water	400 ohm (one face in water)
Typical Motional Resistance Glycerol	3500 ohm (one face in 88 % glycerol)
Sauerbrey Sensitivity Factor	56.6 Hz per ug per cm ² (5 MHz AT-cut at room temperature)
Frequency Mass Coefficient	0.057 Hz per (ng/cm ²) (theory note)
Crystal Holder Material	KynarORing Material: Viton
Holder Connector	BNC
Controller Oscillator Cable	Cat-5, approximately 3 ft
Line Voltage	90 to 132 VAC or 175 to 264 VAC, 47 to 63 Hz
Power Consumption	20 VA total
Operating Mode	analog; external counter and DVM required (not included)
Heavy Load Capability	up to 5000 ohm series resistance
Dry Dial Preset	ten-turn dial 8.0 (dry operation setting per setup procedure)
Varactor Fixed Bias Preset	ten-turn dial 5.0 yields about 6.0 V reverse bias, Cv about 18 pF
Adjust Mode Modulation	75 Hz sine on varactor bias with High/Low/Nullled LEDs
Standard Crystals Included	three chromium/gold QCM crystals

Notes

Specs from SRS QCM100/QCM25 operation manual, QCM100/QCM200 combined catalog, and QCM theory/calibration note. Do not invent counter resolution: frequency accuracy is that of the user-supplied counter and timebase. Forabuilt-in 0.01 Hz counter use QCM200.

Available from Aumictech

Used and refurbished unit evaluated in-house. Calibration: Available on request. Used · In stock.

Request a quote: sales@aumictech.com · +1 (484) 841-9341

Aumictech

788 S 500 E, Pleasant Grove, UT 84062
sales@aumictech.com · +1 (484) 841-9341
<https://aumictechlabs.com>

NIST-traceable calibration · SAM registered ·
30-day returns
NAICS 423690 · 811210 · 541380
Generated August 13, 2026

Technical content from manufacturer specifications. Unit photo/stock from Aumictech listing when available. Specifications subject to OEM documentation and unit condition.